
## MOCON AQUATRAN® 3 - Single Test Report



**Fig. 1:** Water vapor transmission rate of metal oxide barrier on polyimide. The final rate (0.9 mg/m^2-day) can be obtained by subtracting the individual zero (last datapoint) from the last datapoint of the main test (around 310 hours)





Fig. 2: Detail of exposed gold electrode after ALD coating and ultrasonic liftoff of metal oxide barrier.



**Fig. 3:** Accelerated lifetime testing of a preliminary active microelectrode device showing functionality until ~20 days at 87C (~640 days at 37C)



Fig. 4: An ALD-infiltrated device after ALD via opening.